







Legacy Laboratory Services

Outreach Client Swab Guide

Please use the guide below to ensure correct swab use during specimen collection

**The Culture Collect I - Single Swab and Culture Collect II - Double Swab have been discontinued and replaced with the Copan Gene Xpert PCR – Double Swab

Test	Swab
GENITAL	
Chlamydia and Gonorrhoeae by Amplified DNA, BD Qx	E or F
Chlamydia Trachomatis Culture	B
Herpes Simplex Culture	B
Herpes Simplex PCR	B
Streptococcus Group B By PCR	A
Vaginal Pathogens DNA Probe	C
RESPIRATORY	
B Pertussis DNA by PCR	B & D
Beta Strep Rapid Screen	A
Influenza AB RSV PCR	B & D
Nasal Complete MRSA/SA PCR	A
MRSA PCR Nasal Screen	A
Respiratory Panel by PCR	B & D
OTHER	
Bacteria Culture	A
Fungus Culture	A
Culture, Herpes Simplex Virus/Varicella Zoster	A
MRSA Culture	A
MRSA SA Skin PCR Screen	A

Swab	Supply Number
A) Copan Gene Xpert PCR- Double Swab 	274146
B) M-4 Media Viral, Herpes, and Chlamydia Culture 	142081
C) BD Affirm VPIII (Vaginal Pathogens) 	283792
D) Pertussis BBL Culture Swab-Green Top 	248997
E) CT/GC BD Qx Urethral Swab- Male 	308951
F) CT/GC BD Qx Endocervical Swab- Female 	308950