



**LEGACY
LABORATORY
SERVICES**

Pain Medication Compliance Panels (MedManager Screening and Confirmation Summary)

Drug Category	Drug Class	Screening Cutoff	Screening Technique	Confirmation Cutoff	Confirmation Technique
Abused Drugs	Amphetamines	500 ng/mL	EIA	25-50 ng/mL	GC/MS
	Barbiturates	200 ng/mL	EIA	100 ng/mL	GC/MS
	Benzodiazepines	200 ng/mL	EIA	30 ng/mL	LC/MS/MS
	Cocaine	150 ng/mL	EIA	20 ng/mL	GC/MS
	Heroin Metabolite	10 ng/mL	EIA	2 ng/mL	GC/MS
	Marijuana	20 ng/mL	EIA	3 ng/mL	GC/MS
	MDMA	500 ng/mL	EIA	25 ng/mL	GC/MS
	Phencyclidine (PCP)	25 ng/mL	EIA	5 ng/mL	GC/MS
	Alcohol	20 mg/dL	Enzymatic	20 mg/dL	GLC
Opioids	Buprenorphine	20 ng/mL	EIA	NA	NA
	Fentanyl	2 ng/mL	LCMSMS	2 ng/mL	GC/MS
	Methadone	150 ng/mL	EIA	30 ng/mL	LC/MS/MS
	Methadone Metabolite	150 ng/mL	LCMSMS	30 ng/mL	LC/MS/MS
	Opiates	50 ng/mL	EIA	50-100 ng/mL	GC/MS
	Oxycodone	100 ng/mL	EIA	50 ng/mL	GC/MS
	Propoxyphene	150 ng/mL	LCMSMS	30 ng/mL	LC/MS/MS
	Propoxyphene Metab	150 ng/mL	LCMSMS	30 ng/mL	LC/MS/MS
	Tramadol	50 ng/mL	LCMSMS	50 ng/mL	LC/MS/MS
	Tramadol Metab	50 ng/mL	LCMSMS	50 ng/mL	LC/MS/MS
	Meperidine	50 ng/mL	LCMSMS	50 ng/mL	LC/MS/MS
	Meperidine Metab	50 ng/mL	LCMSMS	50 ng/mL	LC/MS/MS

Print new copies from Intranet

Last Approved Date: 04/23/2015

Legacy Laboratory Client Services (503) 413-1234 Toll Free (877) 270-5566