






Legacy Laboratory Services

Outreach Client Stool Collection Guide

Please only provide patients with the appropriate container(s) for test(s) ordered by the provider.

Test	Container
Culture, Stool	A
Gastrointestinal Pathogen Panel by PCR	A
Giardia/Cryptosporidium Detection	B
Parasite Exam Comprehensive	B
Pinworm	C
AFB Culture	D
C. difficile Toxin Detection	D
C. difficile Toxin by PCR	D
Fecal Fat, Qualitative	D
Fecal Leukocytes	D
H. pylori Stool Antigen	D
Norovirus	D
Rotavirus	D
Fecal Immunochemical Occult Blood (FIT) (screening test)	E

<p>A) Enteric Transport Vial</p> 	381820
<p>B) PVA / Formalin Kit</p> 	386129
<p>C) Pinworm Collector Swube Paddle</p> 	220120
<p>D) Sterile 4 oz. Cup</p> 	370167
<p>E) Occult Blood iFOB Collection Mailer Kit</p> 	372280