








Legacy Laboratory Services

Outreach Client Swab Guide

| Test Table Name | Swab |
|-------------------------------------------------------------------------|--------|
| GENITAL | |
| Chlamydia and Gonorrhoeae DNA, BD Qx | F or G |
| Herpes Simplex Virus 1/2 DNA by PCR | D |
| Vaginal Pathogens DNA Probe | H |
| RESPIRATORY | |
| B Pertussis DNA by PCR | D |
| Beta Strep Rapid Screen | B |
| Coronavirus SARS-CoV-2 by PCR | D |
| Influenza AB RSV PCR | D |
| MRSA PCR Nasal Screen | C |
| Nasal Complete MRSA/SA PCR | C |
| Respiratory Panel by PCR | D |
| OTHER | |
| Routine Bacterial Culture (Aerobic Only) | A or B |
| Routine Bacterial Culture (Anaerobic and Aerobic) | E |
| Culture Fungus NOT For Skin, Hair, or Nail Specimens | A or B |
| Culture Methicillin Resistant Staph | A or B |
| MRSA SA Skin PCR Screen | C |

Legacy Laboratory Services provides the supplies listed on this document to our outreach clients. Please refer to our test table for other collection and transport devices that may be acceptable for a given test.

| | | |
|-----------------------------------------|-----------------|--------------------------------------------------------------------------------------|
| A) Culture Collect Single Swab | (Supply 236377) |  |
| B) Culture Collect Double Swab | (Supply 102651) |  |
| C) Transystem Double Swab | (Supply 274146) |  |
| D) M4RT Transport Media | (Supply 142081) |  |
| E) Port-A-Cul Tube | (Supply 102886) |  |
| F) BD ProbeTec Endocervical Swab | (Supply 308950) |  |
| G) BD ProbeTec Urethral Swab | (Supply 308951) | |
| H) BD Affirm VPIII | (Supply 283792) |  |