



Supply Item No. 370167 – Sterile 4 oz Cup

Container	Label Placement Examples
<p>Preferred Patient Label placement:</p> <ul style="list-style-type: none"> ✓ Patient label applied over the manufacturer's barcode label ✓ Applied vertically to container for ease of scanning 	
<p>NOT Preferred Patient Label placement:</p> <ul style="list-style-type: none"> ✗ Too close to manufacturer's barcode which can interfere with scanning the patient's barcode on the label 	
<p>Manufacturer's barcode label</p> <ul style="list-style-type: none"> ✗ Used by the manufacturer to track lots ✗ Do not scan for patient testing 	