



## TEST UPDATES

Notification Date: December 2019

Test Description	Effective Date	Type of Change	Comments
<a href="#">PLATELET P2Y12 TEST</a>	12/3/2019	Deleted	Test deleted
<a href="#">PLATELET ASPIRIN TEST</a>	12/3/2019	Deleted	Test deleted
<a href="#">INFLUENZA AB RSV PCR</a>	12/3/2019	Collect	Updated to: <a href="#">Remel Multi-Microbe Media (M4)</a> is the primary transport media for Legacy Laboratory Services. Legacy will accept these other transport media: <a href="#">Copan Universal Transport Medium (UTM-RT)</a> , <a href="#">BD Universal Viral Transport</a> , <a href="#">Cepheid Xpert Viral Transport media</a> and <a href="#">BD Universal Viral Transport for Quest Diagnostics</a> . See <a href="#">Viral and Pertussis Collection Instructions</a> .
<a href="#">RESPIRATORY PANEL BY PCR</a>	12/3/2019	Collect	Updated to: <a href="#">Remel Multi-Microbe Media (M4)</a> is the primary transport media for Legacy Laboratory Services. Legacy will accept these other transport media: <a href="#">Copan Universal Transport Medium (UTM-RT)</a> , <a href="#">BD Universal Viral Transport</a> , <a href="#">Cepheid Xpert Viral Transport media</a> and <a href="#">BD Universal Viral Transport for Quest Diagnostics</a> . See <a href="#">Viral and Pertussis Collection Instructions</a> .
<a href="#">RHEUMATOID FACTOR</a>	12/11/2019	Handling	Updated to: Allow serum to clot completely at room temperature (minimum: SST - 30 minutes and red top tubes - 60 minutes). Centrifuge and separate serum from cells before transport. Avoid repeated freezing and thawing.
<a href="#">AMYLASE, BLOOD</a>	12/11/2019	Handling	Updated to: Allow serum to clot completely at room temperature (minimum: SST - 30 minutes and red top tubes - 60 minutes). Centrifuge and separate serum from cells before transport.
<a href="#">FLOW CYTOMETRY STUDIES – LEUKEMIA/LYMPHOMA PHENOTYPING</a>	12/11/2019	Name	Updated to: <a href="#">Flow Cytometry Studies – Leukemia/Lymphoma Phenotyping</a>
<a href="#">NICOTINE METABOLITE (COTININE) SEMIQUANTITATIVE SERUM</a>	12/13/2019	Collect	Updated to: Serum or Plasma, one 5.0 mL gold (SST), 8.5 mL red-grey marble, 7.0 mL red, 4.0 mL dark green (sodium heparin), or 4.0 mL lavender (EDTA) top tube

For questions, please call Laboratory Client Services at 503-413-1234 or Laboratory Outreach Client Support at 503-413-4190.



LEGACY  
LABORATORY  
SERVICES

## TEST UPDATES

Notification Date: December 2019

Test Description	Effective Date	Type of Change	Comments
<u>COMPLEMENT C3</u>	12/31/2019	Guidelines  Rejection Criteria	Updated to: <b>Since specimen handling is critical to accuracy of test results, do not order on previously collected specimens (NO "ADD-ONS").</b>  Updated to: <b>Marked icteric, marked hemolyzed, marked lipemic specimens, specimens that have undergone more than one freeze-thaw. No Add-On Requests.</b>
<u>COMPLEMENT C4</u>	12/31/2019	Guidelines  Rejection Criteria	Updated to: <b>Since specimen handling is critical to accuracy of test results, do not order on previously collected specimens (NO "ADD-ONS").</b>  Updated to: <b>Marked icteric, marked hemolyzed, moderate lipemic specimens, specimens that have undergone more than one freeze-thaw. No Add-On Requests.</b>

For questions, please call Laboratory Client Services at 503-413-1234 or Laboratory Outreach Client Support at 503-413-4190.