



**LEGACY  
LABORATORY  
SERVICES**

## TEST UPDATES

**Notification Date: February 2018**

Test Description	Effective Date	Type of Change	Comments
<a href="#"><u>FECAL CALPROTECTIN ELISA</u></a>	2/22/18	Transport, Stability	Stability updated to: Samples may be stored room temperature or refrigerated (2-8°C) for up to 10 days, then frozen at -20°C for 1 year.
<a href="#"><u>UREAPLASMA SPECIES AND MYCOPLASMA HOMINIS CULTURE</u></a>	2/26/18	Mnemonic	Mnemominc was updated to VUREA CULT
<a href="#"><u>UROGENITAL UREAPLASMA AND MYCOPLASMA SPECIES BY PCR</u></a>	2/27/18	New test	
<a href="#"><u>SMEAR, FUNGUS</u></a>	2/27/18	Guidelines, Collect, Handling	
<a href="#"><u>CULTURE, BLOOD FOR FUNGUS</u></a>	2/27/18	Collect, Handling	
<a href="#"><u>INFLUENZA AB RSV PCR</u></a>	2/27/18	Collect	Updated to: Remel Multi-Microbe Media (M4) is the primary transport media for Legacy Laboratory Services. Legacy will accept these other transport media: Copan Universal Transport Medium (UTM-RT, BD Universal Viral Transport, Cepheid Xpert Viral Transport media and BD Universal Viral Transport for Quest Diagnostics.
<a href="#"><u>RESPIRATORY PANEL BY PCR</u></a>	2/27/18	Collect	Updated to: Remel Multi-Microbe Media (M4) is the primary transport media for Legacy Laboratory Services. Legacy will accept these other transport media: Copan Universal Transport Medium (UTM-RT, BD Universal Viral Transport, Cepheid Xpert Viral Transport media and BD Universal Viral Transport for Quest Diagnostics.
<a href="#"><u>URINALYSIS AND MICROSCOPIC WITH CULTURE REFLEX</u></a>	2/27/18	Rejection Criteria	Added: non-sterile containers.
<a href="#"><u>URINALYSIS AND MICROSCOPIC</u></a>	2/27/18	Rejection Criteria	Added: non-sterile containers.



LEGACY  
LABORATORY  
SERVICES

## TEST UPDATES

Notification Date: February 2018

Test Description	Effective Date	Type of Change	Comments
<a href="#"><u>URINALYSIS, MICROSCOPIC ONLY WITH CULTURE REFLEX</u></a>	2/27/18	Rejection Criteria	Added: non-sterile containers.
<a href="#"><u>URINALYSIS, MICROSCOPIC ONLY</u></a>	2/27/18	Rejection Criteria	Added: non-sterile containers.
<a href="#"><u>URINALYSIS, DIPSTICK ONLY</u></a>	2/27/18	Rejection Criteria	Added: non-sterile containers