



LEGACY  
LABORATORY  
SERVICES

## TEST UPDATES

Notification Date: October 2015

Test Description	Effective Date	Type of Change	Comments
<a href="#"><u>ALDOSTERONE/RENIN ACTIVITY RATIO</u></a>	10/12/15	Specimen update; Transport	Deleted heparinized plasma as acceptable specimen type, emphasized specimen must be transported as Critical Frozen
<a href="#"><u>ANTITHROMBIN ACTIVITY</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>C DIFFICILE TOXIN BY PCR</u></a>	10/13/15	Rejection criteria	Formed stools, stool preserved in PVA, SAF or formalin, frozen, specimens past stability
<a href="#"><u>C DIFFICILE TOXIN DETECTION</u></a>	10/28/15	Guidelines, Handling	PCR will be ran for confirmation if either antigen or toxin is positive but not both. Removed "do not freeze" from Handling field.
<a href="#"><u>CHOLESTEROL, LDL, DIRECT</u></a>	10/16/15	Updated Stability	AFTER SEPARATION FROM CELLS: Ambient (18-26°C): 8 hours; Refrigerated (2-8°C): 5 days; Frozen (≤-20°C): 3 months
<a href="#"><u>D-DIMER, QUANTITATIVE, HIGH SENSITIVITY</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>FACTOR VIII: C</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>FECAL CALPROTECTIN ELISA</u></a>	10/27/15	New test	Previously referred to ARUP under test name: CALPROTECTIN, FECAL
<a href="#"><u>FENTANYL SCREEN AND CONFIRMATION, URINE</u></a>	10/13/15	Added Confirmation testing	Confirmatory test will follow if presumptively detected
<a href="#"><u>FIBRINOGEN</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>GAMMA GLUTAMYL TRANSFERASE</u></a>	10/13/15	Added Comments	Legacy Laboratory Services has validated the use of lithium heparin plasma with this assay. Lithium heparin plasma has not been approved or cleared by the U.S. Food and Drug Administration (FDA) for use with this test. Legacy is authorized under Clinical Laboratory Improvement Amendments (CLIA) to perform high-complexity testing.



LEGACY  
LABORATORY  
SERVICES

## TEST UPDATES

Notification Date: October 2015

Test Description	Effective Date	Type of Change	Comments
<a href="#"><u>HEPARIN ASSAY</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>HIV 1/2 ANTIBODY AND P24 ANTIGEN SCREEN WITH REFLEX</u></a>	10/13/15	New test	Replaced: HUMAN IMMUNODEFICIENCY VIRUS 1,2,O COMBINED ANTIBODY SCREEN WITH REFLEX
<a href="#"><u>HYPERCOAGULATION PANEL</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>INFLUENZA AB RSV PCR</u></a>	10/7/15	New test	
<a href="#"><u>IRON BINDING CAPACITY, UNSATURATED (UIBC)</u></a>	10/13/15	Added Comments	Legacy Laboratory Services has validated the use of lithium heparin plasma with this assay. Lithium heparin plasma has not been approved or cleared by the U.S. Food and Drug Administration (FDA) for use with this test. Legacy is authorized under Clinical Laboratory Improvement Amendments (CLIA) to perform high-complexity testing.
<a href="#"><u>LIVER FIBROSIS, FIBRO TEST - ACTITEST PANEL</u></a>	10/1/15	New test	Previously referred to ARUP under test name: HEPATITIS C VIRUS (HCV) FibroSURE
<a href="#"><u>LUPUS ANTICOAGULANT PANEL</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>PAIN MEDICATION COMPLIANCE PANEL, PRESUMPTIVE AUTOCONFIRM</u></a>	10/9/15	New test	
<a href="#"><u>PREALBUMIN</u></a>	10/13/15	Added Comments	Legacy Laboratory Services has validated the use of lithium heparin plasma with this assay. Lithium heparin plasma has not been approved or cleared by the U.S. Food and Drug Administration (FDA) for use with this test. Legacy is authorized under Clinical Laboratory Improvement Amendments (CLIA) to perform high-complexity testing.
<a href="#"><u>PROTEIN C ACTIVITY</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen



LEGACY  
LABORATORY  
SERVICES

## TEST UPDATES

Notification Date: October 2015

Test Description	Effective Date	Type of Change	Comments
<a href="#"><u>PROTHROMBIN TIME-INR</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>PYRUVIC ACID</u></a>	10/13/15	Collection and Stability	Collect Green (sodium or lithium heparin), updated stability to Refrigerated: 2 days
<a href="#"><u>STREPTOCOCCUS GROUP B BY PCR</u></a>	10/28/15	Reported time	Updated to 4 hours from receipt in lab
<a href="#"><u>THROMBIN TIME</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>VON WILLEBRAND PANEL</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>VWF ACTIVITY</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen
<a href="#"><u>VWF ANTIGEN</u></a>	10/28/15	Rejection criteria	Added Whole Blood or Plasma (Refrigerated) as unacceptable specimen