

# **Chlamydia and GC Collection Devices**



# CT/GC BD Qx Endocervical Female Swab

Female Endocervical Specimen Collection Kit for the BD Probe Tec<sup>TM</sup> Chlamydia trachomatis/Neisseria gonorrhoeae (CT/GC) Q<sup>X</sup> Amplified DNA Assays.

### Includes:

Cleaning swab

- White shaft
- Sterile fiber (Rayon) tip

Sterile collection swab

- Scored pink plastic shaft
- Sterile polyurethane foam tip

Swab Diluent for the BD Probe Tec<sup>™</sup> CT/GC Q<sup>X</sup> Amplified DNA Assays

• Black pierceable cap

Storage and handling requirements: 2° to 33°C



# CT/GC BD Qx Urethral Swab

Male Urethral Specimen Collection Kit for the BD Probe Tec<sup>™</sup> Chlamydia trachomatis/Neisseria gonorrhoeae (CT/GC) Q<sup>X</sup> Amplified DNA Assays.

#### Includes:

Sterile collection swab

- Scored blue plastic shaft
- Stainless steel wire
- Sterile rayon tip

Swab Diluent for the BD Probe Tec<sup>TM</sup> CT/GC Q<sup>X</sup> Amplified DNA Assays

Black pierceable cap

Storage and handling requirements: 2° to 33°C



# CT/GC BD Qx Urine Preservative Transport (UPT) tubes

Urine Preservative tubes for the BD Probe Tec<sup>™</sup> Chlamydia trachomatis/Neisseria gonorrhoeae (CT/GC) Q<sup>X</sup> Amplified DNA Assays.

#### Includes:

- UPT tubes
- Transfer pipettes

Storage and handling requirements: 2° to 33°C

Print new copies from Intranet:

Last Approved Date: 05/29/2012

Legacy Laboratory Client Services (503) 413-1234 Toll Free (877) 270-5566