



# SUPPLY REQUEST

Please provide complete information to insure deliverability

**We provide supplies for your use to collect and submit specimens to our laboratories.  
In compliance with Federal regulations, we cannot provide supplies if the specimen is being referred to a non-Legacy laboratory.**

**CLIENT:** \_\_\_\_\_ **SUITE:** \_\_\_\_\_ **CLIENT #:** \_\_\_\_\_ **ORDERED BY:** \_\_\_\_\_

**ADDRESS:** \_\_\_\_\_ **PHONE #:** \_\_\_\_\_

**CITY:** \_\_\_\_\_ **STATE:** \_\_\_\_\_ **ZIP:** \_\_\_\_\_ **DATE:** \_\_\_\_\_

| TUBES  | MICROBIOLOGY COLLECTION   |
|--|---|
| 8.5 ML SST (MARBLE/TIGER TOP) (# OF TUBES) _____ (100/BX) _____ 174130 | BLOOD COLLECTION BOTTLE AEROBIC - SA (EA) _____ 100743            |
| 5 ML SST (GOLD) (# OF TUBES) _____ (100/BX) _____ 101235               | BLOOD COLLECTION BOTTLE ANAEROBIC - SN (EA) _____ 100744          |
| 6 ML RED TOP (PLAIN RED TOP) (# OF TUBES) _____ (100/BX) _____ 101192  | BLOOD COLLECTION BOTTLE PEDIATRIC - PF (EA) _____ 146902          |
| 2 ML LAVENDER (# OF TUBES) _____ (100/BX) _____ 232589                 | BLOOD COLLECTION BOTTLE FUNGUS - FA (EA) _____ 168218             |
| 4 ML LAVENDER (K-EDTA) (# OF TUBES) _____ (100/BX) _____ 158386        | BD SPS SODIUM POLYANETHOLESULFONATE TUBE (EA) _____ 191075        |
| 2.7 ML LTBLUE (COAG) (# OF TUBES) _____ (100/BX) _____ 238095          | FIBER-TIPPED SWAB-STERILE (100/BX) _____ 140762                   |
| 2 ML GRAY (SODIUM FLUORIDE) (# OF TUBES) _____ (100/BX) _____ 342002   | PERTUSSIS BBL CULTURE SWAB - GREEN TOP (EA) _____ 248997          |
| 3 ML PST LITHIUM HEP (# OF TUBES) _____ (100/BX) _____ 147641          | OVA & PARASITE / PARA PAK CULTURE KIT (EA) _____ 204892           |
| 4 ML GREEN (NA-HEPARIN) (# OF TUBES) _____ 135976                      | COMMODE SPECIMEN COLLECTOR (EA) _____ 330514                      |
| 2 ML GREEN (NA-HEPARIN) (# OF TUBES) _____ 158620                      | COPAN GENE XPRT PCR - DOUBLE SWAB (EA) _____ 274146               |
| 6 ML ROYAL BLUE (NO ADDITIVE) (# OF TUBES) _____ 189309                | M-4 MEDIA FOR VIRAL, HERPES & CHLAMYDIA CULTURE (EA) _____ 142081 |
| 6 ML ROYAL BLUE (K2 EDTA) (# OF TUBES) _____ 178793                    | CT/GC BD Qx Endocervical Swab (EA) _____ 308950                   |
| 8.5 ML YELLOW (ACD SOLUTION A) (# OF TUBES) _____ 198556               | CT/GC BD Qx Urethral Swab (EA) _____ 308951                       |
| 6 ML YELLOW (ACD SOLUTION B) (# OF TUBES) _____ 166315                 | CT/GC BD Qx Urine Transport (BX) _____ 308952                     |
| POLYPROPYLENE TRANSFER TUBES 4 ML (EACH) _____ 355643                  | BD AFFIRM VP8 (VAGINAL PATHOGENS) (BX) _____ 283792               |
| CAPS (EACH): RED _____ BLUE _____ LAVENDER _____                       | AMIES ORANGE TOP CHARCOAL SWAB (EA) _____                         |
| YELLOW _____ GREEN _____   | PINWORM COLLECTOR SWAB PADDLE (EA) _____                          |
| QUANTIFERON® GOLD TB COLLECTION KIT (PK) _____ 324781                  | BBL PORT-A-CUL TUBE (EA) _____ 102866                             |

### NEEDLES

|  |
|--|
| 21 GA X 1 1/4" VACUTAINER (SAFETY NEEDLES) (48/BX) _____ 231801  |
| 22 GA X 1 1/4" VACUTAINER (SAFETY NEEDLES) (48/BX) _____ 233726  |
| VACUTAINER HOLDERS (EA) _____ 235763                             |
| SYRINGE SAFETY NEEDLE 23 GA X 1" (305762) (100/BX) _____ 247297  |
| 10 CC SYRINGE LL W/O NEEDLE (50/BX) _____ 247309                 |
| 20 CC SYRINGE W/O NEEDLE (100/BX) _____ 100441                   |
| LATEX-FREE TOURNIQUETS 1" WIDE (EA) _____ 310264                 |
| ALCOHOL WIPES (70% ISOPROPYL) (BX) _____ 100448                  |
| ANTISEPTIC WIPES (NON EDTA) (BX) _____ 310189                    |
| SHARPS CONT. 1.5qt - 8qt. (CIRCLE ONE) (EA) _____ 182079, 218569 |

### MISCELLANEOUS

|   |
|---|
| STERILE 4 OZ CUPS (SEAL W/A CLICK) (100/CS) _____ 135327  |
| 24 HR. COLLECTION BOTTLE (EA) _____ 308880                |
| 24 HR. URINE COLLECTION INSTRUCTIONS (EA) _____ 174565    |
| URINE & STOOL SPECIMEN COLLECTION PAN _____ 356433        |
| GLUCOLA DRINK:  |
| 50 GM LEMON LIME _____ ORANGE _____ (EACH) 145121, 145125 |
| 75 GM LEMON LIME _____ ORANGE _____ (EACH) 224455, 335277 |
| 100 GM ORANGE _____ (EACH) 135503                         |
| HEMOCCULT IFOB COLLECTION MAILER KIT (EA) _____ 306872    |
| VIAL PAP SCI-GEN 3.5 OZ 95% (EA) _____ 197700             |

### CYTOLOGY/PATHOLOGY

|   |
|---|
| SPRAY FIXATIVE, NO 76150 (BOTTLE) _____ 115776                                    |
| CYTOLYTE SOLUTION, NON-GYNECOLOGICAL USE (20/BX) _____ 114359                     |
| SUREPATH CYTOLOGY BROOMS, REMOVAL TIP (25/BAG) _____                              |
| SUREPATH PAP VIALS (COLLECTION DEVICE SEPARATE) (25/PK) _____                     |
| SUREPATH SNAP-OFF BRUSH + SCRAPER (25/PK) _____                                   |
| SUREPATH COMBI BRUSH REMOVAL TIP (25/PK) _____                                    |
| THINPREP CYTOLOGY BRUSHES + SCRAPERS REF# K0012 (25/BAG) _____                    |
| THINPREP PAP VIALS (COLLECTION DEVICE SEPARATE) (20/BX) _____                     |
| THINPREP BROOMS (25/PK) _____   |
| TISSUE BIOPSY BOTTLES (FORMALIN) _____ 40ML 232445 _____ 60ML 150514 (EACH) _____ |
| BBL PORT-A-CUL TUBE (EA) _____ 102866   |

### FORMS / BAGS / LABELS

|  |
|--|
| LABORATORY TEST REQUISITIONS:  |
| _____ CLINICAL LAB # _____   |
| _____ OTHER # _____  |
| SUPPLY REQUEST FORMS (EA) _____ 174515   |
| LOCK BOX: _____ SMALL _____ LARGE _____ BAG (EA) _____ 172972, 230216              |
| SPECIMEN BAGS (DOUBLE POUCH) (100/PKG) _____ 164407                                |
| EUGENE SPECIMEN BAGS (100/PKG) _____ 235707  |
| STICKERS: _____ STAT 341851 _____ PRIORITY 341845 _____ PLASMA 341844 (ROLL) _____ |
| _____ ROOM TEMP 338340 _____ FROZEN 341852 _____ REFRIGERATE 341846 (ROLL) _____   |
| CLEAR PLASTIC BAGS 6 X 9 (100/PK) _____ 143487                                     |
| CLEAR PLASTIC BAGS 12 X 15 (100/PK) _____ 162735                                   |

### ADDITIONAL SUPPLIES

|  |  |
|--|--|
|  |  |
|--|--|

ORDER PICKED BY: \_\_\_\_\_ DATE/TIME: \_\_\_\_\_