

# Legacy Laboratory Services

## Legacy LabAlert

May 2019

An Important Update from Legacy Laboratory Services (LLS)

NEW! Trichomonas Vaginalis by PCR

**TRICHOMONAS VAGINALIS BY PCR**

**Mnemonic: T VAG PCR**

**Supply# 372282, 372283**

LLS now offers Trichomonas Vaginalis testing for both male and female specimens.



Figure A - Supply# 372283

Legacy Laboratory Services has implemented a real-time PCR test for the detection of Trichomonas Vaginalis for the following specimens: **Male and female** urine, endocervical swabs, and vaginal swabs collected under the guidance of health care personnel.




Figure B - Supply # 372282

Specimen Source	Collection Method	Collection Device	Gender
Urine	First-catch Urine	GeneXpert Urine Specimen Collection Kit ( <b>Figure A</b> )	Male and Female
Endocervical/Vaginal	Swab	GeneXpert Vaginal/Endocervical Specimen Collection Kit ( <b>Figure B</b> )	Female

For *technical* information, please contact Megan Draeger at 503-413-2779

For *general* information, please contact your account representative or client services at 877-270-5566



LEGACY  
LABORATORY  
SERVICES

877-270-5566  
[legacylaboratoryservices@lhs.org](mailto:legacylaboratoryservices@lhs.org)